

FREEZER TEMPERATURE LOG – RANGE: $< -15^{\circ}\text{C}$ ($< +5^{\circ}\text{F}$)

Location _____ ECN _____ Month _____ Year _____

*** Daily Sensaphone Checks**

1. The system is plugged in and running. “System ON” light is green.
2. Battery is properly charged or charging. “Battery OK” light is green.
3. The probe is properly positioned inside the unit, in the tube of liquid.
4. Notify Medical Repair if daily checks are not acceptable.

Day of Month	* Daily Sensaphone Check (✓)	AM			PM			Comments
		Time	Temp (°C)	Staff Init	Time	Temp (°C)	Staff Init	
1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								
13								
14								
15								
16								
17								
18								
19								
20								
21								
22								
23								
24								
25								
26								
27								
28								
29								
30								
31								

Corrective Action for out of range temperatures:

1. Store vaccines or other products under proper conditions as soon as possible.
2. Temporarily mark exposed vaccine “do not use”.
3. Call the immunization program (e.g. Milvax) and/or the vaccine manufacturer to determine whether the vaccine is still usable.
4. Document the corrective action on the reverse side of this log.
5. For additional information, or for after hours issues, refer to the flow chart, Enclosure (2), and CDO instructions, Enclosure (3).